hyperSpec Package Review

Devin Gordon

Tuesday, October 30, 2018

#Background

Spectra - What are they?

- "a specific set of values that can vary infinitely within a continuum"
- Generally: 2-dimensional (sometimes higher) dimensional data,
 - where the first dimension identifies points along a continuum,
 - and the second dimension the values of a response of interest
 - at the corresponding points
- Examples: Optical, Mass, Political Alignment

Spectral Data - How do we store them?

- Proprietary file format: not so useful
- .csv files (ASCII): better ~ excel, R
- .spc files (Binary): good ~ R

hyperSpec - How to manipulate .spc files in R

- http://hyperspec.r-forge.r-project.org/
- Handles spectral data
- Features: Convenient import, spectral range selection, shfiting, plotting/viewing, normalization, smoothing, correction, arithmetic, initial data analysis

#hyperSpec Examples

```
##
## Attaching package: 'dplyr'
## The following objects are masked from 'package:stats':
##
##
       filter, lag
## The following objects are masked from 'package:base':
##
##
       intersect, setdiff, setequal, union
## Loading required package: lattice
## Loading required package: grid
## Loading required package: ggplot2
## Package hyperSpec, version 0.99-20180627
##
## To get started, try
      vignette ("hyperspec")
##
##
      package?hyperSpec
      vignette (package = "hyperSpec")
## If you use this package please cite it appropriately.
```


```
## citation("hyperSpec")
## will give you the correct reference.
##
## The project homepage is http://hyperspec.r-forge.r-project.org
##
## Attaching package: 'hyperSpec'
## The following object is masked from 'package:dplyr':
##
## collapse
```

File Import

```
# Read one .spc file and check its structure
# setwd("V:/vuv-data/instr/spectra/")
curfile <- "./data/pet_unstab_cyclic_quv/sa19603_00_ex320.spc"</pre>
step0 <- read.spc(curfile)</pre>
class(step0)
## [1] "hyperSpec"
## attr(,"package")
## [1] "hyperSpec"
step0
## hyperSpec object
      1 spectra
##
##
      4 data columns
      301 data points / spectrum
## wavelength: lambda/nm [numeric] 330.00 331.07 ... 630
## data: (1 rows x 4 columns)
      1. z: x/"a. u." [numeric] 0
##
      2. z.end: x/"a. u." [numeric] 0
##
      3. spc: I/"a. u." [matrix301] 95.5881 98.6830 ... 94.45682
##
      4. filename: filename [character] ./data/pet_unstab_cyclic_quv/sa19603_00_ex320.spc
```

Basic Plotting

```
#Plot with hyperSpec base plotting
plotspc(step0)
```



Multiple .spc Files?

• (Purposefully without a loop or apply)


```
# Read multiple .spc files
\#setwd("./data/pet\_unstab\_cyclic\_quv\_ex320/")
curfile <- "./data/pet_unstab_cyclic_quv/sa19603_00_ex320.spc"</pre>
step0 <- read.spc(curfile)</pre>
curfile <- "./data/pet_unstab_cyclic_quv/sa19603_22_ex320.spc"</pre>
step1 <- read.spc(curfile)</pre>
curfile <- "./data/pet_unstab_cyclic_quv/sa19603_23_ex320.spc"</pre>
step2 <- read.spc(curfile)</pre>
curfile <- "./data/pet_unstab_cyclic_quv/sa19603_24_ex320.spc"</pre>
step3 <- read.spc(curfile)</pre>
curfile <- "./data/pet_unstab_cyclic_quv/sa19603_25_ex320.spc"</pre>
step4 <- read.spc(curfile)</pre>
curfile <- "./data/pet_unstab_cyclic_quv/sa19603_26_ex320.spc"</pre>
step5 <- read.spc(curfile)</pre>
curfile <- "./data/pet_unstab_cyclic_quv/sa19603_27_ex320.spc"</pre>
step6 <- read.spc(curfile)</pre>
curfile <- "./data/pet_unstab_cyclic_quv/sa19603_28_ex320.spc"</pre>
step7 <- read.spc(curfile)</pre>
```

Combine into one object (Good News!)

```
# Use Collpase to combine spectra
spec <- hyperSpec::collapse(step0, step1, step2, step3, step4, step5, step6, step7)</pre>
class(spec)
## [1] "hyperSpec"
## attr(,"package")
## [1] "hyperSpec"
spec
## hyperSpec object
      8 spectra
##
      4 data columns
##
      301 data points / spectrum
##
## wavelength: lambda/nm [numeric] 330.00 331.07 ... 630
## data: (8 rows x 4 columns)
      1. z: x/"a. u." [numeric] 0 0 \dots 0
##
      2. z.end: x/"a. u." [numeric] 0 0 ... 0
##
##
      3. spc: I/"a. u." [matrix301] 95.58810 18.86919 ... 1000
      4. filename: filename [character] ./data/pet_unstab_cyclic_quv/sa19603_00_ex320.spc ./data/pet_un
##
```

Basic Plotting

```
plotspc(spec)
```




Spectral Range Selection

```
#Subset and redifine the hyperSpec object according to wavelength
spec <- spec[,, min ~ 600]</pre>
spec
## hyperSpec object
      8 spectra
##
##
      4 data columns
      271 data points / spectrum
## wavelength: lambda/nm [numeric] 330.00 331.07 ... 600
  data: (8 rows x 4 columns)
      1. z: x/"a. u." [numeric] 0 0 ... 0
##
      2. z.end: x/"a. u." [numeric] 0 0 ... 0
##
##
      3. spc: I/"a. u." [matrix271] 95.58810 18.86919 ... 25.40894
      4. filename: filename [character] ./data/pet_unstab_cyclic_quv/sa19603_00_ex320.spc ./data/pet_un
```

Basic Plotting


```
plotspc(spec)
```



Spectral Arithmetic (Subtraction)


#Perform spectral subtraction

```
sub <- step7 - step0
plotspc(sub[,, min ~ 600])</pre>
```



Removing Bad Data

```
# Let's look at some fake IR data
ir.spc <- chondro</pre>
ir.spc
## hyperSpec object
##
      875 spectra
      5 data columns
##
      300 data points / spectrum
## wavelength: Delta * tilde(nu)/cm^-1 [numeric] 602 606 ... 1798
## data: (875 rows x 5 columns)
##
      1. y: y [numeric] -4.77 -4.77 ... 19.23
      2. x: x [numeric] -11.55 -10.55 ... 22.45
      3. filename: filename [character] rawdata/chondro.txt rawdata/chondro.txt ... rawdata/chondro.txt
##
##
      4. clusters: clusters [factor] matrix matrix ... lacuna + NA
      5. spc: I / a.u. [matrix300] 501.8194 500.4552 ... 169.2942
plotspc(ir.spc, spc.nmax = length(ir.spc))
```



Removing Bad Data


```
#Define any point above 1500 as bad
high.int <- apply(ir.spc > 1500, 1, any)

#Maximum should atleast be 0.1
low.int <- apply(ir.spc, 1, max) < 0.1

#Apply Conditions
ir.spc <- ir.spc[!high.int & !low.int]</pre>
```

Removed the bad spectra

```
plotspc(ir.spc, spc.nmax = length(ir.spc))
```




cq.abs.step0 Correction (One of many options)

```
# Apply basline correction function
blcorr <- spc.fit.poly.below(ir.spc)

# Subtract away correction from original
ir.spc <- ir.spc - blcorr</pre>
```

Corrected Spectra!

```
plotspc(ir.spc)
```



Example of use

Load Cyclic QUV Absorbance Data Set

```
# setwd("v:/vuv-data/proj/3M")
## Read keyfile(s) and set classes

pet.key <- read.csv("./data/pet_unstab_cyclic_quv/3m_sample_key.csv")
pet.key$Sample <- as.character(pet.key$Sample)
pet.key$Product <- as.character(pet.key$Product)
pet.key$Exposure <- as.character(pet.key$Exposure)
pet.key$Step_Retained <- as.character(pet.key$Step_Retained)

## Read step 0 (Unexposed) CyclicQUV Optical Absorbance Data

filenames <- list.files(path = "./data/pet_unstab_cyclic_quv/Cary/CyclicQUV/stepO/", pattern = "\\.spc$

# Read Files

cq.abs.step0 <- lapply(filenames,function(i){
    read.spc(paste("./data/pet_unstab_cyclic_quv/Cary/CyclicQUV/stepO/",i,sep = ""))
})

# Combine resulting list into one hyperSpec object</pre>
```

```
cq.abs.step0 <- hyperSpec::collapse(cq.abs.step0[1:length(cq.abs.step0)])</pre>
# Add step data
cq.abs.step0@data$step <- 0
## Extract sample number from file name using sub (grep) to pull the sample number from the file name m
cq.abs.step0@data$sample <- filenames %>% sub(pattern = "-es00-ms00-mn01_uvvs01.spc",replacement = "")
# Match material info from key file and add to hyperSpec object
rows <- cq.abs.step0@data$sample %>% pmatch(pet.key[,1])
cq.abs.step0@data$material <- pet.key[rows,2]
## Read step 1 CyclicQUV Optical Absorbance Data
filenames <- list.files(path = "./data/pet_unstab_cyclic_quv/Cary/CyclicQUV/step1/", pattern = "\\.SPC$
# Read Files
cq.abs.step1 <- lapply(filenames,function(i){</pre>
  read.spc(paste("./data/pet_unstab_cyclic_quv/Cary/CyclicQUV/step1/",i,sep = ""))
# Combine resulting list into one hyperSpec object
cq.abs.step1 <- hyperSpec::collapse(cq.abs.step1[1:length(cq.abs.step1)])</pre>
# Add step data
cq.abs.step1@data$step <- 1
## Extract sample number from file name using sub (grep) to pull the sample number from the file name m
cq.abs.step1@data$sample <- filenames %>% sub(pattern = "-es01-ms01-mn01_uvvs01.SPC",replacement = "")
# Match material info from key file and add to hyperSpec object
rows <- cq.abs.step1@data$sample %>% pmatch(pet.key[,1])
cq.abs.step1@data$material <- pet.key[rows,2]</pre>
## Read step 2 CyclicQUV Optical Absorbance Data
filenames <- list.files(path = "./data/pet_unstab_cyclic_quv/Cary/CyclicQUV/step2/", pattern = "\\.SPC$
# Read Files
cq.abs.step2 <- lapply(filenames,function(i){</pre>
 read.spc(paste("./data/pet_unstab_cyclic_quv/Cary/CyclicQUV/step2/",i,sep = ""))
})
```

```
# Combine resulting list into one hyperSpec object
cq.abs.step2 <- hyperSpec::collapse(cq.abs.step2[1:length(cq.abs.step2)])</pre>
# Add step data
cq.abs.step2@data$step <- 2
## Extract sample number from file name using sub (grep) to pull the sample number from the file name m
cq.abs.step2@data$sample <- filenames %>% sub(pattern = "-es02-ms02-mn01_uvvs01.SPC",replacement = "")
# Match material info from key file and add to hyperSpec object
rows <- cq.abs.step2@data$sample %>% pmatch(pet.key[,1])
cq.abs.step2@data$material <- pet.key[rows,2]</pre>
Lets check out the contents of each hyperSpec Object
cq.abs.step0
## hyperSpec object
      6 spectra
##
      7 data columns
##
##
      3201 data points / spectrum
## wavelength: lambda/nm [numeric] 1800.0 1799.5 ... 200
## data: (6 rows x 7 columns)
##
      1. z: x/"a. u." [numeric] 0 0 ... 0
      2. z.end: x/"a. u." [numeric] 0 0 ... 0
##
      3. spc: A [matrix3201] 0.1272863 0.1177481 ... 0.2635259
##
      4. filename: filename [character] ./data/pet_unstab_cyclic_quv/Cary/CyclicQUV/step0/sa25000_01-es
##
##
      5. step: [numeric] 0 0 ... 0
##
      6. sample: [character] sa25000.01 sa25002.01 ... sa25010.01
      7. material: [character] B-Melinex243 B-Melinex238 ... B-TeteronU2L92W
cq.abs.step1
## hyperSpec object
##
      54 spectra
##
      7 data columns
      3201 \text{ data points / spectrum}
##
## wavelength: lambda/nm [numeric] 1800.0 1799.5 ... 200
## data: (54 rows x 7 columns)
      1. z: x/"a. u." [numeric] 0 0 ... 0
##
      2. z.end: x/"a. u." [numeric] 0 0 ... 0
##
      3. spc: A [matrix3201] 0.1714494 0.1342367 ... 0.2494305
      4. filename: filename [character] ./data/pet_unstab_cyclic_quv/Cary/CyclicQUV/step1/sa25000_26-es
##
##
      5. step: [numeric] 1 1 ... 1
      6. sample: [character] sa25000.26 sa25000.27 ... sa25010.34
##
      7. material: [character] B-Melinex243 B-Melinex243 ... B-TeteronU2L92W
cq.abs.step2
## hyperSpec object
      48 spectra
```

```
7 data columns
##
##
      3201 data points / spectrum
## wavelength: lambda/nm [numeric] 1800.0 1799.5 ... 200
## data: (48 rows x 7 columns)
      1. z: x/"a. u." [numeric] 0 0 ... 0
##
      2. z.end: x/"a. u." [numeric] 0 0 ... 0
      3. spc: A [matrix3201] -0.030848697 -0.001725956 ... 0.09484347
##
      4. filename: filename [character] ./data/pet_unstab_cyclic_quv/Cary/CyclicQUV/step2/sa25000_27-es
##
##
      5. step: [numeric] 2 2 ... 2
      6. sample: [character] sa25000.27 sa25000.28 ... sa25010.34
##
      7. material: [character] B-Melinex243 B-Melinex243 ... B-TeteronU2L92W
We'll put these all together
# Combine cq.abs.step0s with step 1
cq.abs <- hyperSpec::collapse(cq.abs.step0, cq.abs.step1, cq.abs.step2)</pre>
cq.abs
## hyperSpec object
      108 spectra
##
      7 data columns
      3201 data points / spectrum
## wavelength: lambda/nm [numeric] 1800.0 1799.5 ... 200
## data: (108 rows x 7 columns)
      1. z: x/"a. u." [numeric] 0 0 ... 0
##
      2. z.end: x/"a. u." [numeric] 0 0 ... 0
##
      3. spc: A [matrix3201] 0.1272863 0.1177481 ... 0.09484347
##
      4. filename: filename [character] ./data/pet_unstab_cyclic_quv/Cary/CyclicQUV/step0/sa25000_01-es
##
      5. step: [numeric] 0 0 ... 2
##
      6. sample: [character] sa25000.01 sa25002.01 ... sa25010.34
##
      7. material: [character] B-Melinex243 B-Melinex238 ... B-TeteronU2L92W
```

Utilizing spectral arithmetric to change ordinate axis

```
# Divide by sample thickness for Abs/cm


t.melinex243 <- 0.005
t.melinex238 <- 0.0127
t.melinex626 <- 0.0127
t.melinex618 <- 0.0127
t.teteronu2 <- 0.005
t.teteronu2 <- 0.005
t.teteronu2192w <- 0.0125

for (i in 1:length(cq.abs@data$material)) {

    if (cq.abs@data$material[i] == "B-Melinex243") {cq.abs[i] <- cq.abs[i]/t.melinex243} if (cq.abs@data$material[i] == "B-Melinex238") {cq.abs[i] <- cq.abs[i]/t.melinex238} if (cq.abs@data$material[i] == "A-Melinex626") {cq.abs[i] <- cq.abs[i]/t.melinex626} if (cq.abs@data$material[i] == "A-Melinex618") {cq.abs[i] <- cq.abs[i]/t.melinex618} if (cq.abs@data$material[i] == "B-TeteronU2") {cq.abs[i] <- cq.abs[i]/t.teteronu2} if (cq.abs@data$material[i] == "B-TeteronU2L92W") {cq.abs[i] <- cq.abs[i]/t.teteronu2192W}}</pre>
```

```
cq.abs
## hyperSpec object
##
      108 spectra
##
      7 data columns
##
      3201 data points / spectrum
## wavelength: lambda/nm [numeric] 1800.0 1799.5 ... 200
## data: (108 rows x 7 columns)
      1. z: x/"a. u." [numeric] 0 0 ... 0
##
      2. z.end: x/"a. u." [numeric] 0 0 ... 0
##
      3. spc: A [matrix3201] 25.457266 9.271504 ... 7.587478
##
##
      4. filename: filename [character] ./data/pet_unstab_cyclic_quv/Cary/CyclicQUV/step0/sa25000_01-es
##
      5. step: [numeric] 0 0 ... 2
      6. sample: [character] sa25000.01 sa25002.01 ... sa25010.34
##
      7. material: [character] B-Melinex243 B-Melinex238 ... B-TeteronU2L92W
# Remove errenous points (negative)
`cq.abs` [[`cq.abs` < 0]] <- NA
cq.abs
## hyperSpec object
      108 spectra
##
##
      7 data columns
     3201 data points / spectrum
##
## wavelength: lambda/nm [numeric] 1800.0 1799.5 ... 200
## data: (108 rows x 7 columns)
      1. z: x/"a. u." [numeric] 0 0 ... 0
##
##
      2. z.end: x/"a. u." [numeric] 0 0 ... 0
      3. spc: A [matrix3201] 25.457266 9.271504 ... 7.587478 + NA
##
##
      4. filename: filename [character] ./data/pet_unstab_cyclic_quv/Cary/CyclicQUV/step0/sa25000_01-es
##
      5. step: [numeric] 0 0 ... 2
      6. sample: [character] sa25000.01 sa25002.01 ... sa25010.34
##
     7. material: [character] B-Melinex243 B-Melinex238 ... B-TeteronU2L92W
```

Better plotting using ggplot2 graphics



Peak Information extraction

[1] 133.4993

We can directly address the intensity value at a specific wavelength

```
##
## $y
## [1] 506.0615
fwhm <- upper$y - lower$y</pre>
fwhm
## [1] 122.2706
Imagine this as a function that could operate on multiple spectra
fwhm.spc <- function(spec, peak.wl, lower.bound, upper.bound){</pre>
 peak.wl <- as.numeric(peak.wl)</pre>
   lower <- approx(spec[,,lower.bound ~ peak.wl]$spc, spec[,,lower.bound ~ peak.wl]@wavelength, xout =</pre>
  upper <- approx(spec[,,peak.wl ~ upper.bound]$spc, spec[,,peak.wl ~ upper.bound]@wavelength, xout = c
  fwhm <- upper$y - lower$y</pre>
  center <- mean(c(upper$y,lower$y))</pre>
  skew <- 1 - (upper$y - peak.wl)/(fwhm/2)</pre>
 return(c(fwhm,center,skew))
Let's try it
#Start a data frame to hold the data
ex320data <- read.csv("./data/petpilot-key.csv")
ex320data$Em.Peak <- NA
ex320data$Em.Peak.int <- NA
ex320data$Em.Peak.fwhm <- NA
ex320data$Em.Peak.center <- NA
ex320data$Em.Peak.skew <- NA
ex320data$Em.Peak.skew.mag <- NA
# We'll use a for loop to run each spectrum in the hyperSpec Object
row <- 1
for (i in 1:length(spec)) {
  # Define temporary holder for spectrum of interest
 temp <- spec[i]
  #Define the upper and lower bounds for the approx function
  lower.bound <- min(temp@wavelength)</pre>
    upper.bound <- max(temp@wavelength)</pre>
    # Find the maximum value in the region of interest
    peak <- \max(temp[,,450~470])
    ex320data$Em.Peak.int[row] <- peak
```

```
# Find the peak wavelength

peak.out <- approx(temp[,,450 ~ 470]$spc, temp[,,450 ~ 470]@wavelength, xout = max(temp[,,450~470])
peak.wl <- ifelse(is.na(peak.out$y)==TRUE,450,peak.out$y)
ex320data$Em.Peak[row] <- peak.wl

# Run our function

spec.data <- fwhm.spc(temp, as.numeric(peak.wl), lower.bound, upper.bound)

# Store values in data frame

ex320data$Em.Peak.fwhm[row] <- spec.data[1]
ex320data$Em.Peak.center[row] <- spec.data[2]
ex320data$Em.Peak.skew[row] <- spec.data[3]
row <- row + 1
}</pre>
```

Check the result


```
head(ex320data)
```

```
Material.Type Exposure.Type Exposure.Step
         Sample
## 1 sa19601.00 hydrolytically.stabilized
                                               baseline
## 2 sa19601.01 hydrolytically.stabilized
                                               dampheat
                                                                    1
## 3 sa19601.02 hydrolytically.stabilized
                                               dampheat
                                                                    2
## 4 sa19601.03 hydrolytically.stabilized
                                               dampheat
                                                                    3
## 5 sa19601.04 hydrolytically.stabilized
                                               dampheat
## 6 sa19601.05 hydrolytically.stabilized
                                               dampheat
                                                                    5
    Exposure.Time Total.UV.Dose Em.Peak Em.Peak.int Em.Peak.fwhm
## 1
                              0 450.00
                0
                                           22.99323
## 2
              168
                               0 460.00
                                           266.99850
                                                       122.27058
## 3
                              0 458.93
                                           339.14621
              336
                                                        92.05087
## 4
              504
                              0 460.00
                                           386.45648
                                                         93.64712
## 5
              672
                              0 461.06
                                           408.56882
                                                         90.03078
## 6
              840
                              0 461.96
                                           509.58521
                                                         87.02260
##
    Em.Peak.center Em.Peak.skew Em.Peak.skew.mag
## 1
                NA
## 2
          444.9262 0.246564551
                                               NA
## 3
          458.9813 -0.001115496
                                               NA
## 4
          460.1745 -0.003725775
                                               NA
## 5
          461.9307 -0.019342267
                                               NA
## 6
          463.3233 -0.031332943
                                               NA
```

Other plots


Matrix Plot - Fluorescence Data

```
plotmat(spec)
plotmat(spec, contour = TRUE, add= TRUE)
```




Matrix Plot - IR Data

plotmat(ir.spc)



Stacked Spectra

plotspc(spec, stacked = TRUE)



Summary

hyperSpec

- hyperSpec is a powerful R package for handling and manipulating spectral data in R
- More to come!